

Anti-Nckap1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2A7

Item Number: CM0017213



Introduction

NCKAP1 Rabbit Polyclonal Antibody for WB, ELISA, with cross-reactivity to human, mouse, and rat. Target is part of WAVE complex regulating lamellipodia formation and actin cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	NCKAP1 Rabbit pAb
Remarks/Alias	HEM2; NAP1; NAP125; p125Nap1; NCKAP1
Species	Human
Gene ID (Human)	10787
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human NCKAP1 (NP_038464.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2A7
Protein Weight	129kDa
Shipping	Ice bag