


Anti-Ncoa1 Polyclonal Antibody For Wb, Elisa - Q15788

Item Number: CM0018593



Introduction

High-quality rabbit polyclonal antibody targeting human NCOA1/SRC1, suitable for WB and ELISA. Cross-reacts with mouse. Ideal for studying nuclear receptor coactivation and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	NCOA1 Rabbit pAb
Remarks/Alias	SRC1; KAT13A; RIP160; F-SRC-1; bHLHe42; bHLHe74; NCOA1
Species	Human
Gene ID (Human)	8648
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 720-980 of human NCOA1 (NP_003734.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15788
Protein Weight	157kDa
Shipping	Ice bag