

Anti-Ncoa3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Y6Q9

Item Number: CM0021205



Introduction

NCOA3 Rabbit Polyclonal Antibody (pAb) raised in rabbit against human NCOA3. Validated for WB, IP, ELISA, ideal for studying nuclear receptor signaling, chromatin remodeling, and transcriptional coactivation. - Q9Y6Q9

[Learn More](#)

Parameter	Value
Product Name	NCOA3 Rabbit pAb
Remarks/Alias	ACTR; AIB1; RAC3; SRC3; pCIP; AIB-1; CTG26; SRC-3; CAGH16; KAT13B; TNRC14; TNRC16; TRAM-1; bHLHe42; NCOA3
Species	Human
GeneID (Human)	8202
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NCOA3 (NP_858045.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9Y6Q9
Protein Weight	155kDa
Shipping	Ice bag