


Anti-Ncoa4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13772

Item Number: CM0030442



Introduction

Rabbit polyclonal antibody targeting human NCOA4, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. For ferritinophagy, iron homeostasis, and transcriptional coactivation research.

[Learn More](#)

Parameter	Value
Product Name	NCOA4 Rabbit pAb
Remarks/Alias	RFG, ELE1, PTC3, ARA70, NCOA4
Species	Human
GeneID (human)	8031
GeneID	8031
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NCOA4 (NP_005428.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13772
Protein Weight	70kDa
Shipping	Ice bag