


Anti-Ncoa7 Polyclonal Antibody For Wb, Elisa - Q8Ni08

Item Number: CM0013520



Introduction

High-quality anti-NCOA7 rabbit polyclonal antibody validated for Western blot (WB) and ELISA in human, mouse, and rat samples. NCOA7 is a nuclear receptor coactivator and oxidative stress protector, useful in cancer and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	NCOA7 Rabbit pAb
Remarks/Alias	ESNA1; TLDC4; ERAP140; NCOA7-AS; Nbla00052; Nbla10993; dj187j11.3; NCOA7
Species	Human
GeneID (Human)	135112
GeneID	135112
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human NCOA7 (NP_861447.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NI08
Protein Weight	106kDa
Shipping	Ice bag