


Anti-Ncoa7 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Ni08

Item Number: CM0006094



Introduction

Rabbit monoclonal antibody against human NCOA7, validated for Western blot and ELISA. Recognizes human and mouse NCOA7, a nuclear receptor coactivator involved in transcriptional regulation, oxidative stress protection, and viral entry inhibition.

[Learn More](#)

Parameter	Value
Product Name	NCOA7 Rabbit mAb
Remarks/Alias	NCOA7; ERAP140; ESNA1; NCOA7-AS; Nbla00052; Nbla10993; TLDC4; dj187J11.3; nuclear receptor coactivator 7
Species	Human
Gene ID (Human)	135112
Gene ID	135112
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 777-942 of human NCOA7(NP_861447.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NI08
Protein Weight	12 kDa/25 kDa/28 kDa/33 kDa/92-106 kDa
Shipping	Ice bag