

Anti-Ncor1 Rabbit Polyclonal Antibody For Wb - 075376

Item Number: CM0027680



Introduction

Rabbit polyclonal antibody targeting nuclear receptor corepressor 1 (NCOR1), suitable for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for transcription repression and chromatin dynamics research.

[Learn More](#)

Parameter	Value
Product Name	NCOR1 Rabbit pAb
Remarks/Alias	N-CoR, TRAC1, N-CoR1, hN-CoR, PPP1R109, NCOR1
Species	Human
Gene ID	9611
Immunogen	Recombinant fusion protein containing amino acids 1-200 of human NCOR1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075376
Protein Weight	270kDa
Shipping	Ice bag