


Anti-Ncor2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y618

Item Number: CM0026617



Introduction

Rabbit polyclonal antibody targeting human NCOR2 (SMRT). Suitable for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for studying chromatin remodeling and transcriptional repression mediated by nuclear receptors.

[Learn More](#)

Parameter	Value
Product Name	NCOR2 Rabbit pAb
Remarks/Alias	SMRT; TRAC; CTG26; SMRTE; TRAC1; N-CoR2; TNRC14; TRAC-1; SMAP270; SMRTE-tau; NCOR2
Species	Human
GeneID (Human)	9612
GeneID	9612
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human NCOR2 (NP_006303.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y618
Protein Weight	274kDa
Shipping	Ice bag