

Anti-Ncr1 Rabbit Polyclonal Antibody For Wb And Elisa - O76036

Item Number: CM0028754



Introduction

High-quality rabbit polyclonal antibody against human NCR1 (NKp46, CD335) for Western blot and ELISA. Detects human NCR1; ideal for natural killer cell biology and cancer immunotherapy research.

[Learn More](#)

Parameter	Value
Product Name	NCR1 Rabbit pAb
Remarks/Alias	LY94, CD335, NKP46, NK-p46, NCR1
Species	Human
Gene ID (Human)	9437
Gene ID	9437
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-254 of human NCR1 (NP_004820.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O76036
Protein Weight	34kDa
Shipping	Ice bag