

Anti-Ncr1 Rabbit Polyclonal Antibody For Wb, Elisa - 076036

Item Number: CM0018326



Introduction

Rabbit anti-human NCR1 (CD335, NKp46) polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat natural cytotoxicity triggering receptor 1; ideal for NK cell biology and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	NCR1 Rabbit pAb
Remarks/Alias	LY94; CD335; NKP46; NK-p46; NCR1
Species	Human
GeneID (Human)	9437
GeneID	9437
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 151-303 of human NCR1 (NP_001138929.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76036
Protein Weight	34kDa
Shipping	Ice bag