

Anti-Ncr3 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - O14931

Item Number: CM0007195



Introduction

High-quality Anti-NCR3 Rabbit Monoclonal Antibody validated for WB, IF/ICC, FC (intracellular), and ELISA. Targets natural cytotoxicity triggering receptor 3 (NKp30/CD337) on human NK cells, crucial for tumor cell killing and DC maturation.

[Learn More](#)

Parameter	Value
Product Name	NCR3 Rabbit mAb
Remarks/Alias	NCR3; 1C7; CD337; LY117; MALS; NKp30; natural cytotoxicity triggering receptor 3
Species	Human
GeneID (Human)	259197
GeneID	259197
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC (intra), ELISA
Cross Reactivity	Human
SWISS	O14931
Protein Weight	16kDa/17kDa/18kDa/19kDa/20kDa/21kDa
Shipping	Ice bag