

Anti-Ncs1 Rabbit Polyclonal Antibody For Wb, Elisa - P62166

Item Number: CM0031081



Introduction

Rabbit polyclonal antibody against human NCS1, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for neuronal calcium sensor studies, synaptic plasticity research, and calcium signaling analysis.

[Learn More](#)

Parameter	Value
Product Name	NCS1 Rabbit pAb
Remarks/Alias	FLUP; FREQ; NCS1
Species	Human
GeneID (Human)	23413
GeneID	23413
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human NCS1 (NP_055101.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62166
Protein Weight	22kDa
Shipping	Ice bag