

Anti-Ncs1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62166

Item Number: CM0013382



Introduction

Rabbit polyclonal antibody targeting human NCS1 (Neuronal calcium sensor 1, frequenin homolog). Suitable for WB, IF/ICC, and ELISA applications; cross-reacts with mouse and rat. Predicted molecular weight: 22 kDa.

[Learn More](#)

Parameter	Value
Product Name	NCS1 Rabbit pAb
Remarks/Alias	FLUP, FREQ, NCS1
Species	Human
GeneID (human)	23413
GeneID	23413
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human NCS1 (NP_055101.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62166
Protein Weight	22kDa
Shipping	Ice bag