

Anti-Ncstn Rabbit Pab For Wb, If/Icc, Elisa - Q92542

Item Number: CM0015363



Introduction

Rabbit polyclonal antibody targeting human NCSTN (nicastrin), a gamma-secretase subunit. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for Alzheimer's disease and Notch signaling research.

[Learn More](#)

Parameter	Value
Product Name	NCSTN Rabbit pAb
Remarks/Alias	ATAG1874; NCSTN
Species	Human
Gene ID (Human)	23385
Gene ID	23385
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 429-669 of human NCSTN (NP_056146.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92542
Protein Weight	78kDa
Shipping	Ice bag