


Anti-Ncstn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92542

Item Number: CM0015485



Introduction

Rabbit polyclonal antibody against human NCSTN (Nicastrin, gamma-secretase subunit) for WB, IF/ICC, ELISA. Detects human, mouse, rat. Useful for Alzheimer's disease research and Notch signaling studies.

[Learn More](#)

Parameter	Value
Product Name	NCSTN Rabbit pAb
Remarks/Alias	ATAG1874; NCSTN
Species	Human
GeneID (Human)	23385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-578 of human NCSTN (NP_056146.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92542
Protein Weight	78kDa
Shipping	Ice bag