

Anti-Nctr2 Rabbit Monoclonal Antibody For Wb, Elisa - O95944

Item Number: CM0004614



Introduction

High-quality rabbit monoclonal antibody to human NCTR2 (NKp44, CD336) for WB and ELISA. Reacts with human; ideal for natural killer cell cytotoxicity and tumor immune surveillance research.

[Learn More](#)

Parameter	Value
Product Name	NCTR2 Rabbit mAb
Remarks/Alias	LY95; CD336; NKP44; NK-p44; dj149M18.1; NCTR2
Species	Human
Gene ID (Human)	9436
Gene ID	9436
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NCTR2 (O95944).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95944
Protein Weight	31kDa
Shipping	Ice bag