


Anti-Ndc1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Btx1

Item Number: CM0019712



Introduction

High-quality NDC1 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat NDC1 (~76kDa), a nucleoporin critical for nuclear pore complex assembly. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	NDC1 Rabbit pAb
Remarks/Alias	NET3; TMEM48; NDC1
Species	Human
GeneID	55706
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-530 of human NDC1 (NP_060557.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTX1
Protein Weight	76kDa
Shipping	Ice bag