

Anti-Nde1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nxr1

Item Number: CM0027815



Introduction

Anti-NDE1 Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, rat NDE1 protein. Critical for mitotic spindle orientation and cerebral cortex development.

[Learn More](#)

Parameter	Value
Product Name	NDE1 Rabbit pAb
Remarks/Alias	NDE; LIS4; MHAC; NUDE; NUDE1; HOM-TES-87; NDE1
Species	Human
GeneID (Human)	54820
GeneID	54820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human NDE1 (NP_060138.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXR1
Protein Weight	38kDa
Shipping	Ice bag