

Anti-Ndel1 (Nudel) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Gzm8

Item Number: CM0025464



Introduction

High-quality Anti-NDEL1 (NUDEL) rabbit polyclonal antibody validated for WB and ELISA, with human and mouse cross-reactivity. Ideal for studying microtubule organization, dynein regulation, and neuronal migration. Molecular weight 38kDa.

[Learn More](#)

Parameter	Value
Product Name	NUDEL/NDEL1 Rabbit pAb
Remarks/Alias	EOPA; NDE2; NUDEL; MITAP1; NDE1L1; NUDEL/NDEL1
Species	Human
Gene ID (Human)	81565
Gene ID	81565
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-345 of human NUDEL/NUDEL/NDEL1 (NP_110435.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9GZM8
Protein Weight	38kDa
Shipping	Ice bag