


Anti-Ndp52 Rabbit Monoclonal Antibody For Wb, Elisa - Q13137

Item Number: CM0005544



Introduction

High-quality recombinant rabbit monoclonal antibody against human NDP52 (CALCOCO2), validated for Western blotting and ELISA. Ideal for investigating autophagy and intracellular bacteria degradation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] NDP52 Rabbit mAb
Remarks/Alias	NDP52
Species	Human
GeneID (human)	10241
GeneID	10241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human NDP52 (NP_005822.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13137
Protein Weight	52 kDa
Shipping	Ice bag