


Anti-Ndp52 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13137

Item Number: CM0027700



Introduction

Rabbit polyclonal anti-NDP52 antibody for WB, IF/ICC, ELISA; reacts with human, mouse, rat; detects xenophagy receptor CALCOCO2; molecular weight ~52 kDa; for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	NDP52 Rabbit pAb
Remarks/Alias	NDP52
Species	Human
Gene ID (Human)	10241
Gene ID	10241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human NDP52 (NP_005822.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13137
Protein Weight	52kDa
Shipping	Ice bag