


Anti-NdrG1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92597

Item Number: CM0010294



Introduction

This rabbit monoclonal anti-NDRG1 antibody is validated for Western blot, IHC (paraffin), and ELISA, and cross-reacts with human, mouse, and rat. Ideal for tumor biology and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	NDRG1 Rabbit mAb
Remarks/Alias	GC4; RTP; DRG1; NDR1; NMSL; TDD5; CAP43; CMT4D; DRG-1; HMSNL; RIT42; TARG1; PROXY1; NDRG1
Species	Human
GeneID (Human)	10397
GeneID	10397
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 295-394 of human NDRG1 (Q92597).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92597
Protein Weight	43kDa
Shipping	Ice bag