


Anti-NdrG2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Un36

Item Number: CM0025837



Introduction

Rabbit polyclonal antibody against NDRG2, a tumor suppressor involved in Wnt signaling. Validated for WB and ELISA, reacts with human and mouse. Suitable for cancer and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	NDRG2 Rabbit pAb
Remarks/Alias	SYLD; NDRG2
Species	Human
Gene ID (Human)	57447
Gene ID	57447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 138-357 of human NDRG2 (NP_963832.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UN36
Protein Weight	41kDa
Shipping	Ice bag