


Anti-NdrG3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ugv2

Item Number: CM0013436



Introduction

NDRG3 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human and mouse NDRG3 (41 kDa), a ubiquitously expressed protein with high brain expression. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	NDRG3 Rabbit pAb
Remarks/Alias	NDRG3
Species	Human
GeneID (Human)	57446
GeneID	57446
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 184-363 of human NDRG3 (NP_071922.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UGV2
Protein Weight	41kDa
Shipping	Ice bag