

Anti-NdrG4 Rabbit Polyclonal Antibody For Wb - Q9Ulp0

Item Number: CM0020575

Introduction

Rabbit polyclonal antibody targeting human NDRG4, validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying brain development, ischemic stress, and neuronal signaling.

[Learn More](#)

Parameter	Value
Product Name	NDRG4 Rabbit pAb
Remarks/Alias	BDM1; SMAP8; SMAP-8; NDRG4
Species	Human
GeneID (Human)	65009
GeneID	65009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-220 of human NDRG4 (NP_001229765.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULP0
Protein Weight	38kDa
Shipping	Ice bag