

Anti-Ndst1 Rabbit Polyclonal Antibody For Wb, Elisa - P52848

Item Number: CM0020347



Introduction

High-quality rabbit polyclonal antibody against human NDST1, validated for WB and ELISA. Detects human and mouse NDST1, a key enzyme in heparan sulfate biosynthesis. Ideal for glycobiology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	NDST1 Rabbit pAb
Remarks / Alias	HSST; NST1; MRT46; NDST1
Species	Human
GeneID	3340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-170 of human NDST1 (NP_001534.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52848
Protein Weight	101kDa
Shipping	Ice bag