


Anti-Ndufa1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O15239

Item Number: CM0002622



Introduction

High-quality anti-NDUFA1 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Targets human mitochondrial Complex I subunit, ideal for respiratory chain and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUFA1 Rabbit mAb
Remarks/Alias	MWFE; ZNF183; CI-MWFE; MC1DN12; NDUFA1
Species	Human
GeneID (human)	4694
GeneID	4694
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-70 of human NDUFA1 (O15239).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O15239
Protein Weight	8kDa
Shipping	Ice bag