


Anti-Ndufa1 Rabbit Polyclonal Antibody For Wb, Elisa - O15239

Item Number: CM0027101



Introduction

NDUFA1 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with Human, Mouse, Rat. Targets the mitochondrial Complex I accessory subunit, primarily expressed in heart and skeletal muscle.

[Learn More](#)

Parameter	Value
Product Name	NDUFA1 Rabbit pAb
Remarks/Alias	MWFE; ZNF183; CI-MWFE; MC1DN12; NDUFA1
Species	Human
Gene ID (Human)	4694
Gene ID	4694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human NDUFA1 (NP_004532.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15239
Protein Weight	8kDa
Shipping	Ice bag