


## Anti-Ndufa12 Rabbit Polyclonal Antibody For Wb - Q9Ui09

Item Number: CM0027558



### Introduction

NDUFA12 rabbit polyclonal antibody for WB and ELISA applications, detecting human, mouse, and rat NDUFA12 (17kDa). Ideal for mitochondrial Complex I research and diagnostic assay development.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | NDUFA12 Rabbit pAb                                                                                                 |
| Remarks/Alias    | B17.2; DAP13; MC1DN23; NDUFA12                                                                                     |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 55967                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human NDUFA12 (NP_061326.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q9UI09                                                                                                             |
| Protein Weight   | 17kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |