

Anti-Ndufa4L2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nrx3

Item Number: CM0017208

Introduction

Rabbit polyclonal antibody targeting human NDUFA4L2 (COXFA4L2), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for hypoxia and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	NDUFA4L2 Rabbit pAb
Remarks/Alias	NUOMS; MISTRH; COXFA4L2; NDUFA4L2
Species	Human
GeneID (human)	56901
GeneID	56901
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-87 of human NDUFA4L2 (NP_064527.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRX3
Protein Weight	10kDa
Shipping	Ice bag