


Anti-Ndufa5 Polyclonal Antibody For Wb, Elisa - Q16718

Item Number: CM0025868



Introduction

Rabbit polyclonal antibody against NDUF5, a Complex I accessory subunit. Validated in WB and ELISA for human, mouse, and rat samples. Useful for mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUF5 Rabbit pAb
Remarks/Alias	B13, NUFM, UQOR13, CI-13kB, CI-13KD-B, NDUF5
Species	Human
Gene ID	4698
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-116 of human NDUF5 (NP_004991.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16718
Protein Weight	13kDa
Shipping	Ice bag