

Anti-Ndufa5 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16718

Item Number: CM0015288



Introduction

Anti-NDUFA5 polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects 13kDa subunit of mitochondrial Complex I, highly expressed in heart, muscle, brain.

[Learn More](#)

Parameter	Value
Product Name	NDUFA5 Rabbit pAb
Remarks/Alias	B13; NUFM; UQOR13; CI-13kB; CI-13KD-B; NDUFA5
Species	Human
GeneID	4698 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-116 of human NDUFA5 (NP_004991.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16718
Protein Weight	13kDa
Shipping	Ice bag