


Anti-Ndufa6 Rabbit Polyclonal Antibody For Wb And Elisa - P56556

Item Number: CM0028747



Introduction

High-quality rabbit polyclonal antibody against human NDUFA6, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for mitochondrial complex I research and respiratory chain studies.

[Learn More](#)

Parameter	Value
Product Name	NDUFA6 Rabbit pAb
Remarks/Alias	B14; LYRM6; CI-B14; MC1DN33; NADHB14; NDUFA6
Species	Human
GeneID (Human)	4700
GeneID	4700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-154 of human NDUFA6 (NP_002481.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56556
Protein Weight	15kDa
Shipping	Ice bag