


Anti-Ndufa7 Polyclonal Antibody For Wb, Elisa - O95182

Item Number: CM0024128



Introduction

Rabbit polyclonal antibody against human NDUF7, a mitochondrial complex I subunit. Validated for WB and ELISA, cross-reacts with mouse. Ideal for oxidative phosphorylation and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUF7 Rabbit pAb
Remarks/Alias	B14.5a; CI-B14.5a; NDUF7
Species	Human
GeneID (Human)	4701
GeneID	4701
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human NDUF7 (NP_004992.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O95182
Protein Weight	13kDa
Shipping	Ice bag