

Anti-Ndufa7 Polyclonal Antibody (Ko Validated) For Wb & Elisa - O95182

Item Number: CM0026477



Introduction

Rabbit polyclonal antibody targeting NDUF7 (CI-B14.5a), a mitochondrial Complex I accessory subunit. KO validated for WB and ELISA in human, mouse, and rat samples. For respiratory chain and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NDUF7 Rabbit pAb
Remarks/Alias	B14.5a; CI-B14.5a; A7
Species	Human
Gene ID (Human)	4701
Gene ID	4701
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human NDUF7 (NP_004992.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95182
Protein Weight	13kDa
Shipping	Ice bag