

Anti-Ndufa8 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P51970

Item Number: CM0016938



Introduction

NDUFA8 Rabbit pAb detects the 20kDa mitochondrial Complex I subunit. Validated in WB, IF/ICC, IP, and ELISA with cross-reactivity to human, mouse, and rat. Suitable for respiratory chain and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUFA8 Rabbit pAb
Remarks/Alias	PGIV; CI-19KD; CI-PGIV; MC1DN37; NDUFA8
Species	Human
Gene ID (Human)	4702
Gene ID	4702
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51970
Protein Weight	20kDa
Shipping	Ice bag