


Anti-Ndufa9 Polyclonal Antibody For Wb And Elisa - Q16795

Item Number: CM0025639



Introduction

Rabbit polyclonal antibody targeting NDUF9 (CI-39kD), a mitochondrial Complex I subunit. Validated for WB and ELISA with human, mouse, and rat cross-reactivity. Ideal for mitochondrial and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUF9 Rabbit pAb
Remarks/Alias	CC6; CI39k; COQ11; CI-39k; MC1DN26; NDUFS2L; SDR22E1; NDUF9
Species	Human
GeneID (Human)	4704
GeneID	4704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human NDUF9 (NP_004993.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16795
Protein Weight	43kDa
Shipping	Ice bag