


Anti-Ndufab1 Polyclonal Antibody For Wb, Elisa - O14561

Item Number: CM0017681



Introduction

Rabbit polyclonal antibody against human NDUFAB1, validated for WB and ELISA. Detects endogenous NDUFAB1 in human, mouse, and rat. Suitable for mitochondrial complex I and fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	NDUFAB1 Rabbit pAb
Remarks/Alias	ACP; ACP1; SDAP; FASN2A; NDUFAB1
Species	Human
GeneID (Human)	4706
GeneID	4706
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-156 of human NDUFAB1 (NP_004994.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14561
Protein Weight	17kDa
Shipping	Ice bag