


Anti-Ndufab1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O14561

Item Number: CM0003773



Introduction

NDUFAB1 Rabbit mAb (O14561) validated for WB, IHC-P, ELISA with human, mouse, rat reactivity. Recognizes mitochondrial acyl carrier protein essential in fatty acid synthesis and Complex I; ~17kDa.

[Learn More](#)

Parameter	Value
Product Name	NDUFAB1 Rabbit mAb
Remarks/Alias	ACP; ACP1; SDAP; FASN2A; NDUFAB1
Species	Human
GeneID (human)	4706
GeneID	4706
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-156 of human NDUFAB1 (NP_004994.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14561
Protein Weight	17kDa
Shipping	Ice bag