


Anti-Ndufab1 Rabbit Polyclonal Antibody For If/Icc And Elisa - O14561

Item Number: CM0020783



Introduction

Rabbit polyclonal antibody targeting mitochondrial NDUFAB1, validated for IF/ICC and ELISA, with mouse cross-reactivity. Ideal for investigating Complex I, fatty acid metabolism, and iron-sulfur cluster biogenesis.

[Learn More](#)

Parameter	Value
Product Name	NDUFAB1 Rabbit pAb
Remarks/Alias	ACP; ACP1; SDAP; FASN2A; NDUFAB1
Species	Human
GeneID (Human)	4706
GeneID	4706
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human NDUFAB1 (NP_004994.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	O14561
Protein Weight	17kDa
Shipping	Ice bag