


Anti-Ndufab1 Rabbit Polyclonal Antibody For Wb, Elisa - O14561

Item Number: CM0020898



Introduction

Rabbit polyclonal antibody against human NDUFAB1 (O14561), an acyl carrier protein and complex I subunit. Validated for WB, ELISA; cross-reacts with human, mouse, rat. Ideal for mitochondrial fatty acid synthesis and complex I research.

[Learn More](#)

Parameter	Value
Product Name	NDUFAB1 Rabbit pAb
Remarks/Alias	ACP; ACP1; SDAP; FASN2A; NDUFAB1
Species	Human
GeneID (human)	4706
GeneID	4706
Immunogen	Synthetic peptide: a synthetic peptide corresponding to a sequence within amino acids 50-156 of human NDUFAB1 (NP_004994.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14561
Protein Weight	17kDa
Shipping	Ice bag