

Anti-Ndufaf1 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Y375

Item Number: CM0009659

Introduction

Anti-NDUFAF1 monoclonal antibody validated for WB, IF/ICC, and ELISA. Recombinant rabbit antibody for studying mitochondrial complex I assembly. Reactive with human, 38 kDa. Suitable for mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUFAF1 Rabbit mAb
Remarks/Alias	CGI65; CIA30; CGI-65; MC1DN11; NDUFAF1
Species	Human
Gene ID (Human)	51103
Gene ID	51103
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human NDUFAF1 (Q9Y375)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y375
Protein Weight	38kDa
Shipping	Ice bag