

Anti-Ndufaf2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N183

Item Number: CM0014972

Introduction

Rabbit polyclonal antibody targeting NDUF2 (B17.2L, mimitin), a mitochondrial complex I assembly factor. Validated in WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. For mitochondrial disease and cilia research.

[Learn More](#)

Parameter	Value
Product Name	NDUFAF2 Rabbit pAb
Remarks/Alias	MMTN; B17.2L; MC1DN10; mimitin; NDUFA12L; NDUFAF2
Species	Human
Gene ID (Human)	91942
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-169 of human NDUFAF2 (NP_777549.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N183
Protein Weight	20kDa
Shipping	Ice bag