

Anti-Ndutfaf4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9P032

Item Number: CM0016810

Introduction

Rabbit anti-NDUFAF4 polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for mitochondrial complex I assembly and cancer research.

[Learn More](#)

Parameter	Value
Product Name	NDUFAF4 Rabbit pAb
Remarks/Alias	My013; C6orf66; HRPAP20; HSPC125; MC1DN15; bA22L21.1; NDUFAF4
Species	Human
Gene ID (Human)	29078
Gene ID	29078
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P032
Protein Weight	20kDa
Shipping	Ice bag