


Anti-Ndufb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75438

Item Number: CM0019816



Introduction

NDUFB1 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human and mouse NDUFB1, a subunit of mitochondrial Complex I. Ideal for mitochondrial and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUFB1 Rabbit pAb
Remarks/Alias	MNLL; CI-MNLL; CI-SGDH; NDUFB1
Species	Human
Gene ID (Human)	4707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human NDUFB1 (NP_004536.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75438
Protein Weight	7kDa
Shipping	Ice bag