


Anti-Ndufb10 Polyclonal Antibody For Wb, Ip, Elisa - O96000

Item Number: CM0026626



Introduction

Rabbit polyclonal antibody targeting human NDUFB10, an accessory subunit of mitochondrial complex I. Validated for WB, IP, ELISA in human, mouse, rat. Suitable for respiratory chain and mitochondrial disease studies.

[Learn More](#)

Parameter	Value
Product Name	NDUFB10 Rabbit pAb
Remarks/Alias	PDSW; MC1DN35; NDUFB10
Species	Human
GeneID (Human)	4716
GeneID	4716
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human NDUFB10 (NP_004539.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96000
Protein Weight	21kDa
Shipping	Ice bag