


Anti-Ndufb11 Polyclonal Antibody For Wb, Elisa - Q9Nx14

Item Number: CM0015432



Introduction

NDUFB11 Rabbit pAb is a polyclonal antibody specific for human NDUFB11, a mitochondrial Complex I subunit, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for respiratory chain and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	NDUFB11 Rabbit pAb
Remarks/Alias	ESSS; Np15; P17.3; NP17.3; CI-ESSS; MC1DN30; NDUFB11
Species	Human
Gene ID (Human)	54539
Gene ID	54539
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human NDUFB11 (NP_061929.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NX14
Protein Weight	17kDa
Shipping	Ice bag