


Anti-Ndufb3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43676

Item Number: CM0012872



Introduction

High-quality rabbit polyclonal antibody targeting NDUFB3 (Complex I-B12, 11 kDa), validated for western blot, immunofluorescence/immunocytochemistry, and ELISA. Cross-reacts with mouse. Ideal for exploring mitochondrial respiratory chain function and related disorders.

[Learn More](#)

Parameter	Value
Product Name	NDUFB3 Rabbit pAb
Remarks/Alias	B12; CI-B12; MC1DN25; NDUFB3
Species	Human
Gene ID (Human)	4709
Gene ID	4709
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human NDUFB3 (NP_002482.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	O43676
Protein Weight	11kDa
Shipping	Ice bag