


Anti-Ndufb4 (Ci-B15) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O95168

Item Number: CM0002364



Introduction

NDUFB4 (CI-B15) rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat NDUFB4 protein (~15 kDa). Ideal for mitochondrial complex I research and IVD assay development.

[Learn More](#)

Parameter	Specification
Product Name	NDUFB4 Rabbit mAb
Remarks/Alias	B15; CI-B15; NDUFB4
Species	Human
GeneID (Human)	4710
GeneID	4710
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-129 of human NDUFB4 NP_004538.2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95168
Protein Weight	15kDa
Shipping	Ice bag