


Anti-Ndufb4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95168

Item Number: CM0014661



Introduction

Rabbit polyclonal antibody targeting human NDUFB4, a Complex I subunit. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for mitochondrial and respiratory chain studies.

[Learn More](#)

Parameter	Value
Product Name	NDUFB4 Rabbit pAb
Remarks/Alias	B15; CI-B15; B4
Species	Human
Gene ID (Human)	4710
Gene ID	4710
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human NDUFB4 (NP_001161803.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95168
Protein Weight	15kDa
Shipping	Ice bag